

Product Specification


Product Codes	LX01
Product type	Examination and protective
UDI DI Code	Printed on the box and carton
Chlorinated or Not	Chlorinated
Color	Milky White
Medical grade : (EU) 2017/745	Class I
Protection type: (EU) 2016/425	Category III
Food Contact:	Pass done by SGS
Stretchiness	≥6N
Material	Natural Rubber
Palm thickness:	min. 0.08mm
Finger thickness:	min. 0.08mm
Sterility	Non-Sterile
Texture	Texture at the fingertips
Internal surface	Powdered, Max. 2.0mg/glove
Impermeability level (AQL)	≤1.5
Shape	Ambidextrous
Cuff	Rolled
Glove length	≥240mm
Glove palm width	XS: 70mm, S: 80mm, M: 95mm, L: 110mm, XL: 120mm.
Suitable for sensitive skin	Yes.
Shelf Time:	3 years
storage temperature	5-35°C degrees



Product EAN-13 Codes	XS	S
	 8 859130 638162	 8 859130 638179
M	L	XL
 8 859130 638186	 8 859130 638193	 8 859130 638209

Package Specification

Boxes/Carton:	10
Carton G. Weight:	XS-4.7KG, S-5.1KG, M-5.7KG, L-6.2 KG, XL-7.2 KG
Maximum boxes per EU Pallet:	100
EU Pallets Qty per full truck:	33
EU Pallet G. Weight:	21.5kg
G. Weight per pallet for XS	491.5kg
G.Weight per pallet for S	531.5kg
G.Weight per pallet for M	591.5kg
G. Weight per pallet for L	641.5kg
G. Weight per pallet for XL	741.5kg



Box size: 215mm x 120mm x 80mm

for more information, please visit our website at www.setino.com