

Product Specification

Product Codes	NPF3001P-3005P
Product type	Examination and protective
UDI DI Code	Printed on the box and carton
Chlorinated or Not	Chlorinated
Color	Pink
Protection type: (EU) 2016/425 Category III	
Stretchiness	≥6N
Material	Nitrile
Palm thickness:	min. 0.06mm
Finger thickness:	min. 0.07mm
Sterility	Non-Sterile
Texture	Texture at the fingertips
Internal surface	Powder-Free, Max. 2.0mg/glove
Impermeability level (AQL)	≤1.5
Shape	ambidextrous
Cuff	Rolled
Glove length	≥240mm
Glove palm width	XS: 70mm, S: 80mm, M: 95mm, L: 110mm, XL: 120mm.
Suitable for sensitive skin	Yes. Latex Free
Shelf Time:	5 years
storage temperature	5-35°C degrees



Product EAN-13 Codes	XS	S
	6 948794 412996	6 948794 413009
M	L	XL
6 948794 413016	6 948794 413023	6 948794 413030

Package Specification

Boxes/Carton:	10
Carton G. Weight:	XS-3.8KG S-4.1KG, M-4.4KG, L-4.7KG, XL-5KG
Maximum boxes per EU Pallet:	104 boxes/ pallet
EU Pallets Qty per full truck:	33 pallets
EU Pallet G. Weight:	21.5 kg
G. Weight per pallet for XS	416.7 kg
G.Weight per pallet for S	447.9 kg
G.Weight per pallet for M	479.1 kg
G. Weight per pallet for L	510.3 kg
G. Weight per pallet for XL	441.5 kg



Box size: 220mm x 115mm x 55mm